

FAVARON

CATALOGUE 2015

Equestrian



FAVARON di Favaron Luigi
Via Noventana 123
35027 Noventa Padovana (PD) - ITALY
Tel. +39 049 89 31 288
Fax +39 049 89 31 289
mail favaron@favaron-italy.com
www.favaron-italy.com

*innovative technology***ref. F01 - Front shoe**

Thick: 8mm
 Colour: Neutral
 Material: Rigid

Sizes:
 0—135 x 122mm
 1—148 x 128mm
 2—155 x 135mm
 3—160 x 141mm
 4—169 x 148mm

Remarks:

Good adherence, wear and grip. It keeps the pression uniform on foot.

Sized 0 is rounded at toe for foals shoeing.

ref. F02 - Arrow shoe

Thick: 9mm
 Colour: Neutral
 Material: Rigid

Sizes:
 1—149 x 129mm
 2—150 x 133mm
 3—168 x 149mm

Remarks:

The unique design of the arrow shoe allows the foot to go from landing through weight bearing to breakover with one smooth flowing motion.

ref. F03 - Mushroom shoe

Thick: 10mm
 Colour: Blue
 Material: Rigid / Half-rigid

Sizes:
 1—144 x 120mm
 2—154 x 130mm
 3—167 x 142mm

Remarks:

Take the pressure off the medial and lateral heel and transfer it to the frog: therapeutic device.

ref. F04 - Combi shoe

Thick: 9mm
 Colour: Blue
 Material: Rigid / Half-rigid

Sizes:
 1—149 x 129mm
 2—150 x 133mm
 3—168 x 149mm

Remarks:

The half-rigid version is a perfect device for foals shoeing.

ref. F06 - Spyder shoe

Thick: 3mm
 Colour: Yellow
 Material: Rigid

Sizes:
 0—133 x 120mm
 1—148 x 127mm
 2—155 x 135mm
 3—158 x 140mm
 4—169 x 148mm

Remarks:

Stable, strong and light: it assists the horse in performance to its full potential.

It is like racing without shoes.. Use neckless nails.

ref. FG06 - Spyder with rubber layer

Thick: 5mm
 Colour: Yellow/Black
 Material: Rigid

Sizes:
 0—133 x 120mm
 1—148 x 127mm
 2—155 x 135mm
 3—158 x 140mm
 4—169 x 148mm

Remarks:

Lighter than a composite shoeing, allows a natural biomechanics, protecting the sole like a barefoot horse.

innovative technology**ref. F07 - PG shoe**

Thick: 9mm
 Colour: Neutral
 Material: Rigid

Sizes:
 2—135 x 125mm
 4—150 x 140mm

Remarks:

9mm designed plastic to nail with neckless nails on the horse's foot. Very light: knees saving by absorbing and dispersing ground shocks.


**ref. FG07 - PG with rubber layer**

Thick: 11mm
 Colour: Neutral/Black
 Material: Rigid

Sizes:
 2—135 x 125mm
 4—150 x 140mm

Remarks:

Acts as a shock absorber. Very useful on hard ground, rocks, endurance or street use.

**ref. FINST07 - PG with drilled insert**

Thick: 9mm
 Colours: Yellow for training
 Neutral for racing
 Material: Rigid

Sizes:
 2—135 x 125mm
 4—150 x 140mm

Remarks:

Stable in shape and elastic in movement, a mix of rigidity and flexibility for lasting hoof wellness. plastic mesh allows your horse to wear the bevel at the toe needed.

**ref. F08 - Shoe with studs**

Thick: 8mm
 Colour: Neutral
 Material: Rigid

Sizes:
 2—152 x 135mm

Remarks:

Rigid pins improve grip on hard, iced or mud grounds. Shock absorbing shoes, good wear and durability.

innovative technology

ref. S01 - Flaps




Thick: 9mm
Colour: Neutral
Material: Rigid

Sizes:
1—143 x 122mm
2—150 x 130mm
3—155 x 135mm
4—170 x 150mm

Remarks:
Good traction, rigidity and wear.
Use spyder at toe

ref. S02 - Flaps



Thick: 9mm
Colour: Neutral
Material: Rigid

Sizes:
1—148 x 128mm
2—155 x 135mm
3—160 x 141mm
4—169 x 148mm

Remarks:
Lighter than previous version but easier to use.
Use spyder

ref. S03 - Aired pad




Thick: 4mm
Colour: Neutral
Material: Half-rigid

Sizes:
small —151 x125mm
medium—170 x150mm
large —188 x185mm

Remarks:
It works on the frog providing sole air easing on foot.
Tearing off lateral pressure

ref. S04 - Flat pad thick 3,8mm



Thick: 3,8mm
Colour: Blue
Material: Rigid / Half-rigid

Sizes:
Unique—169 x 145mm

Remarks:
High density plastic protective pad.

ref. S05 - Flat pad thick 2,8mm




Thick: 2,8mm
Colour: Neutral
Material: Soft / Half rigid / Rigid

Sizes:
small —151 x134mm
large —168 x 149mm

Remarks:
Thin, flexible, shock absorbing, resistant.
The elastic material follows heels movement avoiding pad consumption by friction

ref. SG05 - Thin sole with rubber layer



Thick: 4,8mm
Colour: Neutral/Black
Material: Rigid

Sizes:
small —151 x134mm
large —168 x 149mm

Remarks:
Very light shock absorbing

innovative technology




ref. S06 - Flat pad thick 5,5mm

Thick: 5,5mm
Colour: Black
Material: Half rigid

Sizes:
unique—188 x 185mm

Remarks:
Very sturdy, highly tearing resistance for heavy and endurance horses.




ref. S07 - Wedged pad 6/3mm

Thick: 6/3mm
Colour: Blue
Material: Rigid

Sizes:
unique—169 x 145mm

Remarks:
6-3mm graduated for better load distribution




ref. S08 - Wedged pad 3/1mm

Thick: 3/1mm
Colour: Black
Material: Rigid

Sizes:
unique—167 x 145mm

Remarks:
3-1mm graduated for better load distribution



ref. SG08 - Wedged pad 3/1mm with rubber layer

Thick: 5/3mm
Colour: Black
Material: Rigid

Sizes:
unique—167 x 145mm

Remarks:
Therapeutic device



ref. S09 - Sole with magnets

Thick: 3mm
Colour: Neutral
Material:
Rigid / Half rigid

Sizes:
small —151 x134mm
large —168 x 149mm

Remarks:
The new device stimulates blood circulation on foot vascular system. Speeding up recovery time.




ref. S10 - Track

Colour: Neutral
Material: Half rigid

Sizes:
unique—167 x 145mm

Remarks:
Pads with two walls for shoe insert.
Light, ready to use, elastic.
Excellent shock absorbing for gallop horses foot.

innovative technology




ref. S12 - Wall

Colours:
Yellow (rigid)
Green (soft)

Sizes:
small —145 x135mm
large —165 x 145mm

Remarks:
The internal hole allows better grip on ground.




ref. S13 - Wedged sole 10/3mm

Thick: 10/3mm
Colour: Blue
Material: Rigid

Sizes:
unique—165 X145mm

Remarks:
It modifies foot's positioning aiding pumping blood.




ref. S14 - Twice

Thick: 5,2mm
Colour: Yellow
Material: Rigid

Sizes:
unique—167 x 143mm

Remarks:
Rigid sole working on heels tearing off pressure on frog.
Helpful for horses needing more ground coverage on foot back bottom.



ref. S15 - Bull

Thick: 6mm
Colour: Yellow
Material: Rigid

Sizes:
unique—165 x 145mm

Remarks:
Soft sole used to give frog protection or frog support, stiffen a shoe that has tendency to spread.
Reduce the excessive load and trauma on crushed or compromised heels enhancing healing.

innovative technology

PRODUCTS SELECTION

Ref. F07 - PG shoe

Designed plastic to nail with neckless nails on the horse's foot.
Very light: knees saving by absorbing and dispersing ground shocks.



Features

Thick: 9mm
Colour: Neutral
Material: Rigid

Sizes:
2—135 x 125mm
4—150 x 140mm

Ref. FG07 - PG with rubber layer

Acts as a shock absorber.
Very useful on hard ground, rocks, endurance or street use.



Features

Thick: 11mm
Colour: Neutral/Black
Material: Rigid

Sizes:
2—135 x 125mm
4—150 x 140mm

Ref. FINST07 - PG with drilled insert

Stable in shape and elastic in movement, a mix of rigidity and flexibility for lasting hoof wellness.
plastic mesh allows your horse to wear the bevel at the toe needed.



Features

Thick: 9mm
Colours:
Yellow for training
Neutral for racing
Material: Rigid

Sizes:
2—135 x 125mm
4—150 x 140mm

Ref. F06 - Spyder shoe

Stable, strong and light: it assists the horse in performing to its full potential.
It is like racing without shoes.
Rubber version is lighter than a composite shoeing: it allows a natural biomechanics, protecting the sole like a barefoot horse.
Use neckless nails.



Features

Thick: 3mm.

Colour: Yellow/Black
Material: Rigid

Rubber version
Thick: 5mm.

Sizes:
0—133 x 120mm
1—148 x 127mm
2—155 x 135mm
3—158 x 140mm
4—169 x 148mm

HELP

Have a nice and safe walk with
HELP (with no metal parts)

**ECONOMIC, DURABLE,
LIGHT AND COMPACT**

The boot "Help" offers an alternative for horses health treatments topical remedies, poultices and pre-medicated pads and shoe loss. The one size "Help" is designed for fitting the majority of horses. The fastening of the boot "Help" ensures a safe walk and it is metal free.

The rockered urethane outsole provides comfort, traction and avoids slipping on all surfaces.

The boot "Help" can be inserted one in the other in order to avoid dirty to get in the boot. It is easy to remove , to clean and to put on again

**A PRECIOUS HELP IN THE
THERAPEUTIC TREATMENTS**



art. Help

Help is suitable for sizes from 4 to 7



Closed version



Version with central hole

Open version